

XTERN CLASSIC+



The XTERN Classic+ is
geTTing an upgrade



overhauled design

THE XTERN SUMMIT'S GENIUS SEE-THROUGH DESIGN HAS BEEN ADAPTED TO THE XTERN CLASSIC+.



new "drill-less" design

PUT DOWN THE DRILL! THE XTERN CLASSIC+ CAN NOW BE ASSEMBLED WITHOUT ANY POWER TOOLS OR ADDITIONAL HOLES, THANKS TO ITS NEW MONOCORE TECHNOLOGY® IT NOW BOASTS A ONE-PIECE BODY AND NEW CALF SECTION.



upgraded calf cuff

THE CALF CUFF IS NOW HEIGHT ADJUSTABLE AND EQUIPPED WITH PRESSURE SCREWS FOR EASIER ASSEMBLY.



updated all-in-one instruction manual

OUR NEW, REDESIGNED ASSEMBLY GUIDE IS NOW MORE USER FRIENDLY WITH PICTOGRAMS AND MINIMAL TEXT.

+91 7835880155

DROP FOOT MADE SIMPLE

MAXIMIZING FUNCTIONALITY & COMFORT



The TurboMed XTERN CLASSIC+ AFO is a revolutionary new treatment option for drop foot. It's unique patented design acts as an exoskeleton to the impaired limb and provides the user with unparalleled levels of function without discomfort or rubbing.

The XTERN CLASSIC+ is a custom build modular ankle foot orthosis (AFO) made from highly durable thermoplastic (limited warranty of 2 years).

It attaches to the outside of footwear and is quickly and easily interchangeable from shoe to shoe.

Unlike many other AFO's, the XTERN does not prevent ankle plantar flexion or limit dorsiflexion, making it easier and much more comfortable.

FEATURES & BENEFITS:

- Exoskeleton Design Fits onto the Outside of the Shoe
- No Skin Contact with Foot or Ankle
- Energy Return due to patented design
- Easily Conforms to Uneven Ground and Slopes
- Highly Durable Thermoplastic (2 year warranty)
- Heat Adjustable (to change width and dorsiflexion strength)
- Easily Transferable from Shoe to Shoe
- Symmetrical Design, Fits both Right & Left Foot
- Easily Positioned to Invert or Evert the Foot if Required

CLINICAL INDICATIONS:

- Drop Foot & Dorsiflexion Weakness
- Peroneal Nerve Injury
- Charcot Marie Tooth (CMT)
- Stroke (CVA)
- Multiple Sclerosis (MS)
- Cerebral Palsy (CP)
- Guillain-Barre Syndrome (GBS)
- Motor Neuropathy
- Motor Neurone Disease (MND)

CAN BE USED WITH:

- Trainers/Running Shoes
- Lace and Velcro-Fastening Shoes/Boots
- Walking Shoes/Boots
- Safety Work Boots/shoes
- Dress Shoes
- Sandals

CLINICAL LIMITATION

- Inappropriate Footwear
- Gross Fixed Deformities of the Foot and Ankle
- Severe Spasticity (Hypertonicity)
- Quadriceps weakness
- Severe Proximal weakness
- Crouch gait



Turbo  **Med**[™]
ORTHOTICS

OUR WEBSITE

